

## Bur 系列

(直機)



### A. ROUND BURS

Item Number	Head Shank	Head Material	Head Diam.	Bur Number	Total Length	Cutting Length	Features
21-44-0060	HP	Carbide	1.8mm	6	44mm	-	10/Pkg.
21-50-0060	HP	Carbide	1.8mm	6	50mm	-	10/Pkg.
21-65-0060	HP	Carbide	1.8mm	6	65mm	-	10/Pkg.
21-44-0080	HP	Carbide	2.3mm	8	44mm	-	10/Pkg.
21-50-0080	HP	Carbide	2.3mm	8	50mm	-	10/Pkg.
21-65-0080	HP	Carbide	2.3mm	8	65mm	-	10/Pkg.

(直機)



### B. CYLINDER BUR

40-50-0016	HP	Carbide	1.5mm	-	50mm	7.5mm	10/Pkg.
------------	----	---------	-------	---	------	-------	---------

(直機)



### C. TAPER BURS(CARBIDE)

61-44-0002	HP	Carbide	1.6mm	702	44mm	5mm	10/Pkg.
61-50-0002	HP	Carbide	1.6mm	702	50mm	5mm	10/Pkg.
85-21-4886	HP	Carbide	1.2mm	701	44mm	5mm	10/Pkg.
85-21-4883	HP	Carbide	1mm	700	44mm	5mm	10/Pkg.

(直機)



### D. OVAL BONE CUTTING BURS

31-68-0031	HP	Carbide	4mm	-	68mm	9mm	Cross Cut
31-50-0023	HP	Carbide	5.5mm	-	50mm	10mm	Cross Cut
504-080505	HP	Carbide	5.5mm	-	70mm	10mm	Cross Cut

(直機)



### E. LINDEMANN BURS

19-26-0162	FG	SS	1.6mm	-	26mm	9mm	-
10-44-0166	HP	SS	2mm	-	44mm	10mm	-
10-44-0166S	SN	SS	2mm	-	44mm	10mm	-
10-65-0168	HP	SS	2.3mm	-	65mm	24mm	-
10-65-0168S	SN	SS	2.3mm	-	65mm	24mm	-
10-80-0169	HP	SS	2.3mm	-	80mm	35mm	-

(彎機)




### F. KIRSCHNER BURS (前端有鋸齒狀可修正角度)

88-33-0010	CA	SS	3.15-4mm	-	34mm	15mm	Internal Irr.
88-33-0009	CA	SS	1.6-2.35mm	-	34mm	15mm	Internal Irr.

**HOT**  
精選

(彎機)



### G. ROUND BURS FOR BONE CUTTING (含內沖修骨Bur)

11-26-0020	CA	SS	2mm	-	26mm		Int.Irrigation
11-27-0015	CA	SS	4mm	-	34mm		Int.Irrigation
11-27-0017	CA	SS	6mm	-	34mm		Int.Irrigation
V-1-811C030RAXL	CA	SS	3mm	-	34mm		Int.Irrigation

**HOT**  
精選



### Q 6966 Sterilization basket

made of stainless steel, for small parts,  
直徑 46mm

(彎機)



### H. CRESTOTOME BUR (修整骨平面Bur)

11-26-0357	CA	SS	3.75-7mm	-	26mm	7mm	Internal Irr.
------------	----	----	----------	---	------	-----	---------------



1

### I. Sinus Scoring Diamond Burs

直機、彎機手機專用



 編號：W-9025-0018

 編號：W-90440018

2

3

### J. Sinus Lift Egg Burs

45mm, 2.8HD

 編號：W-00450028

4

5

### K. Stop Pilot Twist Drill

小直徑的Drill

Ø 1.2mm, Single Patient-Sterile

 MDLSPD

Ø 1.2mm, Single Patient 13mm Long-Sterile

 MDLSPDL

Ø 1.5mm, Stop 8-10-11, 5-13-15mm

 D1515 D1515

Ø 2.0mm, Stop 8-10-11, 5-13-15mm

 D2015 D2015

6

7

8

9

10

### L. IntraLock™ Drill

三角定位

Lancer 0-2mm-Sterile

 編號：LD

### M. Step Drill (Lag Screw)

三角定位



 編號：V-186RF018-1

### N. Step Drill (Lag Screw)

補骨打洞專用

0.8mm, 1.6mm Steps, Contra-Angle

 W-454-281001

Ø 1.1mm, Optional Pilot Drill, CA


 W-454-40201

Ø 0.7mm, Bur Twist 25mm  
6mm Stop Contra-Angle


 W-502-50610



#### Stop Pilot Twist Drills-Ø1.5mm




#### Stop Pilot Twist Drills-Ø2.0mm





## ACE 取骨環鑽





### Trephine Burs - Contra Angle

一般長度 (彎機專用)

Manufactured from titanium (Ti-6Al-4V ELI).

Item number	Shank	Head Material	I.D.	O.D.	Total Length	Cuting Length	Markings
A 04-9480-01	CA	Ti	2mm	2.8mm	35mm	18.5mm	8,10,13,15,18mm
B 04-9480-02	CA	Ti	2.7mm	3.5mm	35mm	18.5mm	8,10,13,15,18mm
C 04-9482-01	CA	Ti	3.75mm	4.6mm	35mm	18.5mm	8,10,13,15,18mm
D 04-9482-02	CA	Ti	4.25mm	5.10mm	35mm	18.5mm	8,10,13,15,18mm
E 04-9484-01	CA	Ti	5.25mm	6.10mm	32mm	15.5mm	8,10,13,15mm
F 04-9485-01	CA	Ti	6mm	6.9mm	32mm	15.7mm	8,10,13,15mm
G 04-9485-02	CA	Ti	7mm	7.9mm	32mm	15.7mm	8,10,13,15mm
H 04-9487-01	CA	Ti	8mm	8.9mm	32mm	14mm	8,10,13mm
I 04-9487-02	CA	Ti	10mm	10.95mm	32mm	14mm	8,10,13mm



### DR. Smiler's Trephine Burs - Posterior

後牙區專用較短淺

Manufactured from titanium (Ti-6Al-4V ELI).



Item number	Shank	Head Material	I.D.	O.D.	Total Length	Cuting Length	Markings
J 04-9486-01	CA	Ti	6mm	6.9mm	23mm	6.75mm	6mm
K 04-9486-02	CA	Ti	7mm	7.9mm	23mm	6.75mm	6mm



### Latch To Square Driver Converter

L 20-1554-11 14mm long, secures and locks contra shaft for hand or ratchet drive


彎機改手動轉接頭



### Trephine Burs - Handpiece

直機專用

Manufactured from titanium (Ti-6Al-4V ELI).



Item number	Shank	Head Material	I.D.	O.D.	Total Length	Cuting Length	Markings
M 04-9492-01	HP	Ti	3.75mm	4.6mm	52mm	18.5mm	8,10,13,15,18mm
N 04-9492-02	HP	Ti	4.25mm	5.1mm	52mm	18.5mm	8,10,13,15,18mm
O 04-9495-01	HP	Ti	6mm	6.9mm	49mm	15.7mm	8,10,13,15mm
P 04-9497-01	HP	Ti	8mm	8.9mm	49mm	14mm	8,10,13mm
Q 04-9497-02	HP	Ti	10mm	10.95mm	49mm	14mm	8,10,13mm


#### 臨床操作圖





### Kohler™ Trephine Burs

彎、直機取骨環鑽 材質：不銹鋼 stainless



彎機手機使用




編號	內徑	外徑
Q-6942	2.5mm	3.0mm
Q-6521	3.0mm	3.5mm
Q-6522	5.0mm	6.0mm
Q-6523	7.5mm	8.5mm
Q-6524	10mm	11mm

直機手機使用



編號	內徑	外徑
Q-6742	4.0mm	5.0mm
Q-6743	6.0mm	7.0mm
Q-6744	10mm	11mm



Q-7109 長度3.2cm, 內徑 Ø 3.5mm

### Kohler™ Rotary tissue punch (for contra-angles)

彎機手機軟組織切割器



編號	Q-7297	Q-7298	Q-7299	Q-7300	Q-7301	Q-7316	Q-7317
外徑	3.0mm	3.5mm	4.1mm	4.5mm	5.1mm	6.1mm	7.1mm
CMU							



編號	RPCA3	Q-7356	Q-7357	Q-7358
外徑	3.0mm	4.0mm	5.0mm	6.0mm

### Intralock™ Tissue Push

彎機手機軟組織切割器




### Kohler™ Rotary tissue punch

彎機手機軟組織切割器 (前端有定位點)

with centering pin (for contra-angles)



彎機使用



## Kohler™ Circular knife 手動式定位軟組織切割器

by hand, with centering pin



定位



外徑 MUC	3.0mm	3.5mm	4.1mm	4.5mm	5.1mm	6.1mm	7.1mm
	Q-7291	Q-7292	Q-7293	Q-7294	Q-7295	Q-7296	Q-7302

## Kohler™ Tissue punches KEYES 手動軟組織切割器

Tissue punches KEYES (直)



Q-7516	長度8 cm, 內徑 Ø 3mm
Q-7517	長度8 cm, 內徑 Ø 4mm
Q-7518	長度8 cm, 內徑 Ø 5mm
Q-7519	長度8 cm, 內徑 Ø 6mm
Q-6427	長度10 cm, 內徑 Ø 6mm
Q-7520	長度8 cm, 內徑 Ø 8mm

Tissue punches KEYES, 90度 (彎)



Q-7521	長度8 cm, 內徑 Ø 3mm
Q-7522	長度8 cm, 內徑 Ø 4mm
Q-7523	長度8 cm, 內徑 Ø 5mm
Q-7524	長度8 cm, 內徑 Ø 6mm
Q-6874	長度10 cm, 內徑 Ø 6mm
Q-7525	長度8 cm, 內徑 Ø 8mm

## Kohler™ Tissue punches KEYES, MINI 手動軟組織切割器

Tissue punches KEYES, MINI



Q-7500	長度3.5cm, 內徑 Ø 3.5mm
Q-7502	長度3.5cm, 內徑 Ø 5.0mm
Q-7504	長度3.5cm, 內徑 Ø 6.5mm
Q-7506	長度3.5cm, 內徑 Ø 8.0mm

Tissue punches KEYES, MINI 1/2



Q-7501	長度3.5cm, 內徑 Ø 3.5mm
Q-7503	長度3.5cm, 內徑 Ø 5.0mm
Q-7505	長度3.5cm, 內徑 Ø 6.5mm
Q-7507	長度3.5cm, 內徑 Ø 8.0mm



Q 6966 Sterilization basket

made of stainless steel, for small parts,  
直徑 46mm













## Disposable Tissue punches

拋棄式軟組織切割器 Sterile, Single Use, Plastic Handle

W-603-3331A	Ø 1.5mm	W-603-3334	Ø 4mm
W-603-3331B	Ø 2.5mm	W-603-3335	Ø 5mm
W-603-3332	Ø 3mm	W-603-3336	Ø 6mm
W-603-3333	Ø 3.5mm		




直徑對照表1:1

1.5mm	2mm	2.5mm	3mm	3.5mm	4mm
					
5mm	6mm	7mm	8mm	9mm	10mm
					






# ACE Surgical Blades & Handles 軟組織切割器



For Taking Precise Keratinized Epithelial And Connective Tissue Grafts

### W-65-000-14 Unigraft Knife Handle Only

The same handle is used for the ACE Unigraft Knife and the ACE Connective Tissue Graft Harvester. Replacement Shoes are Interchangeable.



取皮方向 →

Step(1)

Step(2)

Step(3)



作一深切2-3mm深度之V形溝槽，以方便取皮刀取皮。  
★深切之深度請依取皮刀之規格而定

## ●可依您想要的取皮厚度選擇刀片，使臨床更加順手



### Single Blade Replacement Shoe 單刀片

· Ideal for Epithelial Grafts

編號	取皮厚度 X 寬度
W-65-000-04	1.25mm x 6mm
W-65-000-05	1.5mm x 6mm
W-65-000-07	1.25mm x 9mm
W-65-000-08	1.5mm x 9mm








### Double Blade Replacement Shoe 雙刀片

· Ideal for Connective Tissue and Epithelial Grafts

編號	取皮厚度 X 寬度
W-65-000-09	1.5mm x 9mm

## TIPS FOR SUCCESSFUL USE OF THE CONNECTIVE TISSUE GRAFT HARVESTER



**Graft size 1.5 by 6mm or 9mm depending on the blade width used.**

The harvested piece of connective tissue has a consistent 1.5mm thickness and a 9mm width. The length is controlled by the ACE CT Graft Harvester advance and is 10mm for this example.


特色：

1. 單刀片取出為單層表皮
2. 雙刀片手術即可切出 2 片

### 取皮臨床操作圖



## Kohler™ Scalpel blades and holders 刀柄系列



一次可裝2片刀片  
兩刀片間距 1.0mm 兩刀片間距 1.5mm

180°


依所需角度隨時調整  
可靈活運用於臨床中  
★萬向刀柄  
附17cm & 13cm  
長、短刀柄各一支

Q-4209 16 cm  
Q-4210 16 cm  
Q-5523 3  
Q-6442 with 2 handles 13/17cm  
Q-5527 16 cm for 2 blades  
Q-5528 16 cm for 2 blades

Q-5527, Q-5528操作方式

## Handle for micro scalpel blades

進口手術專用刀片 規格：100片/盒



Micro scalpel blades

Q-5032  
Q-5034

Q-6747  
Handle for micro scalpel blades  
Ø 6.5 mm

Q-7693  
Ø 6.5 mm

#11 一般刀片 #12 一般刀片 #15 一般刀片

◎ 比一般刀片直徑較小

#12d雙面刀片 W-0036212  
#15C細頭微彎 W-0037315