



整組工具明細如下：

- 1 W-454-19001 truFIX Driver Handle 握柄
- 2 W-454-17001 truFIX Driver Tip 圖釘專用-手動接頭
- 3 W-454-15001 truSCREW Driver Tip 螺絲專用-手動接頭
- 4 W-454-15201 truSCREW Driver Tip, Contra-Angle 螺絲專用-彎機接頭
- 5 W-454-17901 truSCREW Driver Removal Sleeve
- 6 W-50-27-00045 Pilot Bur, 0.45mm, 27mm Length, Contra-Angle
- 7 W-454-40201 1.1mm Pilot Drill, Optional, CA Drill 螺絲專用鑽針
- 8 W-454-281001 Step Drill(Lag Screw), 0.8, 1.6, 2.1mm Steps, Contra-Angle
- 9 W-454-50201 truSCREW Packaging Removal Tool
- 10 W-20-1577-02 0.88mm External Hex Hand Driver
- W-454-18501 truFIX Tray 工具收納盒



truSCREW™ truTACK™ 規格表

truSCREW™ 1.2 Thread Diameter

| | | | | |
|-------------|----------|-----|----------------|-------|
| W-454-27202 | truSCREW | 3HD | 1.2TD x 4.5TL | 5/pkg |
| W-454-27203 | truSCREW | 3HD | 1.2TD x 6.0TL | 5/pkg |
| W-454-27204 | truSCREW | 3HD | 1.2TD x 7.5TL | 5/pkg |
| W-454-27205 | truSCREW | 3HD | 1.2TD x 9.0TL | 5/pkg |
| W-454-27206 | truSCREW | 3HD | 1.2TD x 10.5TL | 5/pkg |

truSCREW™ 1.5 Thread Diameter

| | | | | |
|-------------|----------|-----|----------------|-------|
| W-454-27502 | truSCREW | 3HD | 1.5TD x 6.0TL | 5/pkg |
| W-454-27503 | truSCREW | 3HD | 1.5TD x 7.5TL | 5/pkg |
| W-454-27504 | truSCREW | 3HD | 1.5TD x 9.0TL | 5/pkg |
| W-454-27505 | truSCREW | 3HD | 1.5TD x 10.5TL | 5/pkg |
| W-454-27506 | truSCREW | 3HD | 1.5TD x 12.0TL | 5/pkg |
| W-454-27507 | truSCREW | 3HD | 1.5TD x 13.5TL | 5/pkg |
| W-454-27508 | truSCREW | 3HD | 1.5TD x 15.0TL | 5/pkg |

truSCREW™ 2.0 Thread Diameter


| | | | | |
|-------------|----------|-----|----------------|-------|
| W-454-28002 | truSCREW | 3HD | 2.0TD x 6.0TL | 5/pkg |
| W-454-28003 | truSCREW | 3HD | 2.0TD x 7.5TL | 5/pkg |
| W-454-28004 | truSCREW | 3HD | 2.0TD x 9.0TL | 5/pkg |
| W-454-28005 | truSCREW | 3HD | 2.0TD x 10.5TL | 5/pkg |
| W-454-28006 | truSCREW | 3HD | 2.0TD x 12.0TL | 5/pkg |
| W-454-28007 | truSCREW | 3HD | 2.0TD x 13.5TL | 5/pkg |
| W-454-28008 | truSCREW | 3HD | 2.0TD x 15.0TL | 5/pkg |

TENTING SCREW 1.5 Thread Diameter

| | | | | |
|-------------|--------------|-----|----------------|--------|
| W-454-25501 | TentingScrew | 4HD | 1.5TD x 7.0TL | 1 each |
| W-454-25502 | TentingScrew | 4HD | 1.5TD x 11.0TL | 1 each |
| W-454-25503 | TentingScrew | 4HD | 1.5TD x 15.0TL | 1 each |



truTACK™

| | | | | |
|-------------|---------|-------|--------|--------|
| W-96-003-13 | truTACK | 2.5HD | 3mm TL | 10/pkg |
| W-96-003-14 | truTACK | 2.5HD | 5mm TL | 10/pkg |



獨特八角設計不易掉落，臨床使用更順手

truSCREW™ 2.5HD truTACK™ SCREW





1

2

3

4

5


6

7


8

9


10









Insert the bone screw



Screw the Self-Drilling Surgical Bone Screw into the recipient site.



TITANIUM MESH GRAFTING WITH ACE BONE SCREWS




註：此病例流程以可吸收再生膜為例



使用ACE TruFix BoneScrew固定bone block



補上 Bio-oss、FDDBA、Unigraft 骨粉



不同選擇的 Tenting Screw






FlatOne Tenting Screw



FlatOne tenting screw 工具規格表

| 產品代碼 | 產品規格 |
|-----------|--|
| FABS4 | FlatOne bone screw $\phi 2 \times 4\text{mm}$ |
| FABS5 | FlatOne bone screw $\phi 2 \times 5.5\text{mm}$ |
| FABS7 | FlatOne bone screw $\phi 2 \times 7\text{mm}$ |
| FABS8 | FlatOne bone screw $\phi 2 \times 8.5\text{mm}$ |
| FABS10 | FlatOne bone screw $\phi 2 \times 10\text{mm}$ |
| FABS11 | FlatOne bone screw $\phi 2 \times 11.5\text{mm}$ |
| TIFANPCS1 | FlatOne cover screw (H : 1.2mm) |
| Y-001 | Hex Driver for FlatOne cover screw-0.038inch |
| IFACAD | FlatOne bone screw driver (C. A) |
| DL16 | FlatOne bone screw hand driver (long) |
| DS16 | FlatOne bone screw hand driver (short) |






ACE Bone Screw Tacking Instruments 鈦合金二合一工具組(舊系統套件表)



舊包裝



454-19501 Small Handle, Titanium, 3-1/2", 9cm 握柄-小



454-19001 Large Handle, Titanium, 4", 10cm 握柄-大

1.0mm Spline 1mm接頭適用於2.5mm HD(直徑)



454-15501 1mm Spline Driver Tip, 55mm 手動專用接頭-短



454-15502 1mm Spline Driver Tip, 85mm 手動專用接頭-長



454-15701 1mm Spline Driver Tip, CA, 19mm 接彎機手機專用接頭

1.25mm Spline 1.25mm接頭適用於3.0-4.0mm HD(直徑)




454-15001 1.25mm Spline Driver Tip, 55mm 手動專用接頭-短



454-15002 1.25mm Spline Driver Tip, 85mm 手動專用接頭-長



454-15201 1.25mm Spline Driver Tip, CA, 19mm 接彎機手機專用接頭



W-96-003-15 Tack Seating Instrument, Titanium, Straight, Includes Cover 握柄(直)



W-96-003-16 Tack Seating Instrument, Titanium, Angled, Includes Cover 握柄(有角度)



W-64-003-12S Surgical Mallet, Stainless Steel 錘子

Available Screw Sizes

自鑽/自攻式骨固定螺絲規格 材質：鈦合金



Micro Screws - 1.2mm thread diameter - 包裝：5pc/set
2.5mm Head Diameter (Use the 1.0mm Spline Driver Tip)

- 454-21201 1.2mm thread diameter length 3.0mm
- 454-21203 1.2mm thread diameter length 6.0mm



Micro Screws - 1.5mm thread diameter - 包裝：5pc/set
2.5mm Head Diameter (Use the 1.0mm Spline Driver Tip)

- 454-21304 1.5mm thread diameter length 7.5mm
- 454-21305 1.5mm thread diameter length 9.0mm



Mid Screws - 1.5mm thread diameter - 包裝：5pc/set
3.0mm Head Diameter (Use the 1.25mm Spline Driver Tip)

- 454-21501 1.5mm thread diameter length 4.5mm
- 454-21502 1.5mm thread diameter length 6.0mm
- 454-21507 1.5mm thread diameter length 13.5mm
- 454-21508 1.5mm thread diameter length 15.0mm



Mini Screws - 2.0mm thread diameter - 包裝：5pc/set
3.0mm Head Diameter (Use the 1.25mm Spline Driver Tip)

- 454-22002 2.0mm thread diameter length 7.5mm
- 454-22003 2.0mm thread diameter length 9.0mm
- 454-22004 2.0mm thread diameter length 10.5mm
- 454-22005 2.0mm thread diameter length 12.0mm



Drill Point Tenting Screws 單支包裝
4.0mm Head Diameter (Use the 1.25mm Spline Driver Tip)

- 454-26502 2.3mm thread diameter length 11.0mm

Step Pilot Drills 鑽針



454-44001 Step Pilot Drill, CA, SS, 30mm, (Step Diameter:0.8,1.6,2.1mm)



454-44501 Step Pilot Drill, HP, SS, 48mm, (Step Diameter:0.8,1.6,2.1mm)



W-50-44-00045 Twist Drill, 0.45mm diam., 4mm CL, HP, 44mm TL 鑽針(直機)



W-50-27-00045 Twist Drill, 0.45mm diam., 4mm CL, CA, 27mm TL 鑽針(彎機)



W-20-1577-01 0.88mm External Hex Screwdriver, 73mm (Not Included In The Kit) 圖釘專用起子




ACE Micro Osseo-Plate System 鈦骨板系列


產品比例為1:1

Micro Osseo-Plates 使用1.2mm與1.5mm的Bone Screws

Titanium



ACE Titanium Self-Drilling Screws for Micro Osseo-Plate System (0.5mm High)




配合Bone Screw使用，詳細介紹請參考P2-6

A27 頭顱模型



編號：5-U-3B-A27 十分解
編號：5-U-3B-A20 一般標準三分解





Ti-Mesh 美國原裝進口

衛署醫器輸字第001400號

Ti-MESH 鈦網特殊型



- A** W-003-9434
7 x 14mm, 0.2mm thick



Ti-MESH Plate 鈦網圓孔型



Mesh Plate

- C** W-003-9430 120 x 60mm, 0.1mm thick
- D** W-003-9433 152 x 66mm, 0.2mm thick



Ti-MESH 鈦網 (小)

- B** W-003-9429
34 x 25mm, 0.1mm thick(單片包裝)



- Fully resorbable
- Easy to handle
- Does not adhere to itself when hydrated
- Zero memory
- 3-6 month predictable resorption

conFORM™
Resorbable Collagen MEMBRANE

ZERO MEMORY

Memorable Handling.



1

2

3

4

5

6

7

8

9

10

8





骨再生的最佳選擇利器

CYTOFLEX® (金色-鈦網)

衛署醫器輸字第017451號

Ti-Mesh 0.8mm Thick Precision Titanium Mesh



1 C02-4401
25 x 30 mm




2 C02-4501
30 x 40 mm



3 C02-4301
20 x 25 mm



4 C02-4201
12 x 25 mm



Small Membrane Punch 打洞專用再生膜鑷子

可打洞可定位



W-96-003-03 3", 7.5cm, Stainless Steel

Membrane Holding Forceps 夾取再生膜專用鑷子

Q-6862
W-96-003-22 Without Lock, Angled



W-96-003-24 With Lock


Membrane Placement Instrument



ACE Buser Type Elevator for Membrane Placement and Periosteal Elevation
W-03-001-39 7", 18cm



Q-3167 Raspatories, 24G, 17.5cm



Q-3171 BUSER, 17.5cm



CYTOFLEX®

衛署醫器輸字第017451號

銷售十年以上 · 美國原裝進口

www.unicarebiomedical.com

TEF-Guard™ Membrane 微孔隙設計

Microporos PTFE Membrane


不可吸收 再生膜

類似Gore-Tex Membrane

微隙/多孔性設計，可提供血液最佳供應來源，並提供軟組織附著防止再生膜易位之優點，並有效阻隔細菌侵入。



編號：C04-0105 (凹面)
規格：12 x 24mm (5片/盒)



編號：C04-0205 (凹面)
規格：25 x 30mm (5片/盒)



編號：C01-0105 (平面)
規格：25 x 30mm (5片/盒)

臨床操作圖



術前PANO



植入IC5011STI Implant



放骨粉及不可吸收再生膜 (TEF-GUARD)



覆蓋再生膜



傷口無法完全縫合



一星期後拆線



手術後3個星期抽再生膜



再生膜特寫




再生膜抽出傷口特寫



抽出再生膜一星期後拆線



術後追蹤-傷口已癒合



術後PANO

1

2

3

4

5

6

7

8

9

10

GBR GTR & 骨再生的最佳選擇利器

CYTOFLEX® 含鈦支架再生膜

衛署醫器輸字第017451號


銷售十年以上 · 美國原裝進口




Ti-Enforced

Ti Reinforced Membrane
含鈦支架再生膜

NEW 第四代




1 C05-0101
13 x 20mm




2 C05-0301
11 x 21mm




3 C05-0501
14 x 24mm




4 C05-0701
18 x 29mm




5 C05-0901
19 x 26mm




6 C05-1101
27 x 29mm



7 C05-1501
36 x 25mm




8 C05-1901
41 x 30mm



9 C05-1301
20 x 25mm




10 C05-1701
25 x 30mm



11 C05-2101
32 x 40mm



先利用包裝紙測量補骨區大小
再剪裁 Ti-reinforce membrane





在台銷售10年以上 品質值得信賴

UNIGRAFT™ Synthetic Bioactive glass

美國原裝進口

衛署醫器輸字第016661號

合成人造骨

合成材料，全素患者可適用。
最佳造骨骨粉，可與任何材料共同使用，提供骨細胞最佳養分。



適用於：

1. 植牙最佳補骨材料
2. 牙周病骨缺陷修補
3. 拔牙後的骨缺陷修補
4. 可與任何骨粉混合使用，達到養分供給
5. 使血小板凝結達到骨細胞快速生長

編號：P-E-U03-0105 規格：500 mg / 瓶 1000 mg / 瓶
 P-E-U03-0205 Conveniently packaged in 0.5 gram vials,
 P-E-U03-0305 5 per box, easy to use.

下顎前牙補骨 臨床病例



Geistlich Bio-Oss 吉士德 百歐士

衛署醫器輸字第020576號



Small granules(0.25-1mm)
 細顆粒(S)
 Quantities : 0.5g



Small granules(1-2mm)
 粗顆粒(L)
 Quantities : 0.5g

GBR GTR & 骨再生的最佳選擇利器

ReOss® SYNTHETIC NAND-HA INFUSED BIOMATERIAL

奈米結構人工合成生物補骨材料

衛部醫器輸字第026322號


ReOss® 奈米結構 | 6-12個月 HA Particles 生長 新生骨



最佳使用：

HA + TCP 材料

1. 口內部分骨缺損
2. 牙周骨間或骨缺損
3. 齒槽脊缺損、拔牙窩
4. 上顎竇增高手術




ReOss® Powder

規格：RP-1 0.5 cc
RP-3 1 cc




ReOss® Putty

規格：RP-2 0.5 cc
RP-4 1 cc



放大5,000倍電顯圖，可以看到 ReOss 具多孔性的PLGA基質。



放大50,000倍電顯圖，可以看到 ReOss 的 HA顆粒吸附在PLGA基質上。

注入新活力的奈米化 人工合成鈣磷灰石 生物材料

ReOss Putty + Bio-Oss





一般外科專用 大面積止血材料

Johnson & Johnson 快速止血紗布(可吸收性)

衛署醫器輸字第004219號

快速止血紗布可吸收性 (Surgicel)

特色：

- 1.快速止血
- 2.用於取皮臨床手術
- 3.下顎取骨塊可填傷口

編號：P-C-B1903GB (24片/盒)



SURGICEL止血氧化再生纖維：

是由pure alpha cellulose所製造的止血棉，能在微血管、小靜脈、小動脈及其他滲血的部位促進血小板的附著並形成凝塊以達到止血的效果。

其特性：

- 2至8分鐘達到止血（另有"Nu-Knit"厚片規格，只須1至5分鐘即可達成止血）
- 7至14天內完全被人體吸收。
- 含殺菌性-對Gram Positive / Negative微生物有殺菌性。
- 以植物為原料，可減少動物或人體病源的污染。



規格：5cm x 7cm



門診拔牙 最佳使用止血紗布

可溶性止血癒合紗布

衛署醫器製字第003479號

特色：

- 1.拔牙後傷口癒合使血凝塊，可填滿傷口達到快速骨生長
- 2.撐鼻竇時可用小規格搭配使用

- 3至5分鐘達到止血，5-10日內完全被人體吸收
- 止血是它的功能，幫助傷口癒合才是它的效果
- 快速止血、快速吸收、快速癒合
- 高分子聚合物高黏度，體表結合性強不易脫落
- 吸收組織液提供快速癒合環境
- 塊狀、圓狀包裝，不需剪裁，符合您所需的規格

編號P-C-B-C-001：圓狀裝 0.5 x 0.7CM 60PCS/盒

編號P-C-B-C-002：塊狀裝 1.5 x 1.5CM 60PCS/盒



編號：P-C-B-C-002



編號：P-C-B-C-001

Bone Wax 骨臘

填補空間材料

Made From Beeswax And Othe Natural Materials

衛署醫器輸字第006172號

編號：P-D-AW810



2.5克/個 (24個/盒)

GBR & GTR

骨再生的最佳選擇利器

拔牙後加速傷口癒合 · 輔助植牙的最佳利器 衛署醫器製字第003043號



第4代

規格：10mm x 20mm

編號：SA001020

速吉安® SurgiAid® Collagen Wound Dressing 膠原蛋白傷口敷料



- ★ 吸收期 6-8 週
- ★ 牛筋製成 (澳洲牛)
- ★ 柔軟可壓縮
吸水率高 恢復性快
- ★ 夾取容易 手術操作方便
- ★ 可濕式或乾式使用
- ★ 連通孔洞結構 支撐强度高

優越的臨床操作性 · 高強度的拉伸特性

衛署醫器製字第003195號

敷美安® FormaAid® Collagen Membrane 可吸收再生膜

- 吸收期長達16-18週
- 牛筋製成 (澳洲牛)
- 高抗拉強度及抗穿刺
- 可帶針縫合
- 臨床操作方便
- 良好的服貼性



實體拍攝
高抗拉強度



第4代

| 編號 | 規格 |
|----------|---------------------|
| FA-T3040 | 30mm x 40mm x 0.3mm |
| FA-T2530 | 25mm x 30mm x 0.3mm |
| FA-T1520 | 15mm x 20mm x 0.3mm |



1

2

3

4

5

6

7

8

9

10




美國原裝進口

ACE Micro Saws

All ACE Micro Saws feature:

- Slim saw blades keep the loss of bone to a minimum
- Low noise level and high cutting precision
- Wide field of application
- Use with any E-Coupling electric motor, 1:1 ratio
- Actual Size: 5-1/2", 14cm
- Applications: Maxillofacial Surgery, Implant Dentistry, Craniofacial Surgery, Plastic Surgery, Hand Surgery, Podiatry
- Includes Fork Wrench



E-COUPLING: Use with any E-Coupling electric motor, 1:1 ratio, with speeds 15-20,000 RPM

W-682402


ACE Micro Sagittal Saw

- Can make horizontal or vertical bone cuts
- Does not cut soft tissue – cuts bone aggressively
- As the blades are extremely thin (0.35mm), very delicate cuts can be achieved
- Ideal for ridge splitting or block grafts
- Recommended Maximum Speed: 20,000 RPM

ACE Micro Sagittal Saw Blades

- Blade Thickness 0.35mm
- Blades are extremely thin, very delicate cuts can be achieved

| | 編號 | 角度 | 寬度 | 長度 |
|----------|-----------|------------------|------|------|
| A | W-682470 | Straight | 4mm | 10mm |
| B | W-682470A | 45° Angle | 4mm | 10mm |
| C | W-682471 | Straight | 6mm | 10mm |
| D | W-682471A | 45° Angle | 6mm | 10mm |
| E | W-682475 | Straight, Round | 6mm | 10mm |
| F | W-682475A | 45° Angle, Round | 6mm | 10mm |
| G | W-682472 | Straight | 6mm | 15mm |
| H | W-682472A | 45° Angle | 6mm | 15mm |
| I | W-682473 | Straight | 10mm | 15mm |
| J | W-682474 | Straight | 12mm | 27mm |







W-682400
ACE Micro Oscillating Saw

• Recommended Maximum Speed: 15,000 RPM

ACE Micro Oscillating Saw Blades

- Blade Thickness 0.3mm
- Makes Vertical or Horizontal Cuts
- Saw Blades with Short Shanks, Working Length is 11mm
- Saw Blades with Long Shanks, Working Length is 26mm

| 尺寸 | 編號 | 長度 |
|----------------------|----------|-------------|
| A 6mm x 14mm | W-682424 | Short, 27mm |
| | W-682427 | Long, 42mm |
| B 10mm x 14mm | W-682425 | Short, 27mm |
| | W-682428 | Long, 42mm |
| C 15mm x 14mm | W-682426 | Short, 27mm |
| | W-682429 | Long, 42mm |
| D 6mm x 9mm | W-682430 | Short, 27mm |
| | W-682431 | Long, 42mm |




W-682405
ACE Micro OsseoScalpel

• Recommended Maximum Speed: 20,000 RPM

ACE Micro OsseoScalpel Saw Blades

Blade Thickness 0.35mm

| 編號 | 長度 |
|----------|------|
| W-682435 | 12mm |
| W-682436 | 18mm |
| W-682437 | 24mm |



W-682410
ACE Micro Compass Saw/Power Periosteal



This Micro Compass Saw can be used as a Power Periosteal simply by changing the blade!

- Recommended Maximum Speed as a Micro Compass Saw: 15,000 RPM
- Recommended Speed as a Power Periosteal: 500-5,000 RPM Maximum

ACE Micro Compass Saw Blades

Blade Thickness 0.35mm

| 編號 | 長度 |
|----------|------|
| W-682412 | 11mm |
| W-682413 | 14mm |
| W-682414 | 18mm |
| W-682415 | 22mm |
| W-682416 | 26mm |





1

2

3

4

5

6

7

8

9

10



彎機/手動 兩用擴骨組


ACE Ridge Expanding Kit

編號：W-0556-0000 (整組編號)



整組工具明細：

- 1 W-055-60201 2.5mm Bone Expander
- 2 W-055-60301 3.1mm Bone Expander
- 3 W-055-60401 3.6mm Bone Expander
- 4 W-055-60501 4.3mm Bone Expander
- 5 W-04-9734-01 Connection Driver Tool (Hex) for CA Handpiece
- 6 W-96-002-28 Universal Ratchet
- 7 W-20-1558-11 Driver Extender (Hex)
- 8 W-96-004-49 Square Titanium Bowl
- 9 W-20-1950-11 Autoclavable Tray 收納盒





ACE RIDGE EXPANDING CHISEL KIT

劈骨專用器械組




Markings at 8,10,13,15,18mm





W-454-4700 Set of 4 Ridge Chisels 整組編號

整組工具明細：

- W-454-4701 3.0mm Ridge Chisel
 - W-454-4801 4.0mm Ridge Chisel
 - W-454-4901 5.0mm Ridge Chisel
 - W-454-4501 9.5mm Ridge Chisel
 - W-765-1090 Autoclavable Tray 收納盒
 - W-454-1501 Chisel Handle (Uses 1.25mm Hex Driver) 蝴蝶扳手
- All chisels and tray may be purchased separately.




W-454-1501 蝴蝶扳手
Chisel Handle Optional, With Attachment Screw
(Use 1.25mm Hex Driver)





Bone Harvesting Instrument



ACE Sagittal Split Osteotomy Chisels

取下顎骨專用側彎 Chisels

Useful in elevating the proximal segment



側面微彎

W-0502017 Cottle Nasal Chisel 6mm, 6-1/4", 16cm 微彎



W-0501018 8mm Chisel, Flat, 45° Angle, 6", 15cm

W-0501220 5mm Chisel, Flat, 45° Angle, 6", 15cm

Spatula Osteotomy Instruments 取骨塊專用 Chisels



W-0400069 4mm x 1mm, Slight Curve, 8", 20cm



W-0500010 4mm x 30mm, 6", 長15cm 超薄 (厚度約0.45mm)



W-0500011 5mm x 35mm, 6", 長15cm 超薄 (厚度約0.75mm)

1

2

3

4

5

6

7

8

9

10



Periodontal Knives 牙周專用刮匙



3-HU-KGF8
8
Goldman-Fox

3-HU-KGF11
11
Goldman-Fox

3-HU-KB5/6
W-417-9480
5/6
Buck

3-HU-KS1/2
W-417KS12
1/2
Sanders

Probe



3-6-9-12

1-2-3-4-5-6-7-8-
9-10-11-12-13-
14-15

12-9-6-3

3-HU-PQ2N
Q-2N
Nabers
Color coded

CP-15UNC
University of
North Carolina


另備有Hu-Friedy相關器械目錄，歡迎來電洽詢



1



Periodontal Knives 牙周專用刮匙



3-HU-KK15/16
W-417-9468
15K/16K
Kirkland

3-HU-KGF7
7
Goldman-Fox

3-HU-KSOLT1/2
1/2
Solt

3-HU-K01/2
1/2
Orban

3-HU-SPR1/2
PR1/2
Prichard
DE Curette

Suture Pliers




3-HU-SP20
20
CORN
Suture Pliers
15.2cm/6"

另備有Hu-Friedy相關器械目錄，歡迎來電洽詢



Periodontal Files

牙周專用刮匙



3-HU-FS9/10S
9/10
Schluger Curved File

3-HU-FB11/12
11/12
Buck Straight File


3-HU-FS1/2S
1S/2S
Sugarman File M/D

3-HU-FS3/4S
W-417-9439
3S/4S
Sugarman File B/L

3-HU-S13K/TG
13K/TG

3-HU-S13K/13KL
13K/13KL

另備有Hu-Friedy相關器械目錄，歡迎來電洽詢



Periodontal Chisels 牙周專用刮匙

Prichard

Buser



3-HU-CKN1/2
1/2
Kramer-Nevins

3-HU-C36/37
36/37
Rhodes Back Action

3-HU-S44/45
44/45
Remington Hoe

3-HU-PPR3
PR-3
Prichard
Periosteal

3-HU-PPBUSER
Buser

3-HU-CRM
Retractor U
Of Minnesota

另備有Hu-Friedy相關器械目錄，歡迎來電洽詢



Periotome Blades 牙周病專用剝離器

立即拔牙植牙手術專用器械好幫手！



3-HU-CIGFT4
#4
Extra-Flex

3-HU-CIGFT3
#3
Extra-Flex



Periodontal Instruments 植牙周圍炎專用刮匙 (可替換)



3-HU-IMPL4R/4L
Columbia

3-HU-IMPLH6/H7
Plasteel Tip/Handle
(一組有配10個)



3-HU-IMPL4R/4L(12個)
W-417-4403(10個)
4R/4L
Plasteel Tips 塑膠材質



3-HU-IMPL204S(12個)
W-I204S(10個)
I204S
Plasteel Tips 塑膠材質



3-HU-IMPL6/7(12個)
W-IH6/7(10個)
H6/H7
Plasteel Tips 塑膠材質



3-HU-IMPHDL6
W-417-4406
double-end
#6
Satin Steel

另備有Hu-Friedy相關器械目錄，歡迎來電洽詢




1



Gracey Curette

Mini Five



3-HU-SAS1/2
1/2 Gracey

3-HU-SG11/12
11/12

3-HU-SG13/14
13/14

3-HU-SG15/16
15/16

3-HU-SG3/4
3/4

3-HU-SG7/8
7/8

磨刀石



3-HU-SS3C
No.3
Ceramic
Medium grit ^{w/d}
7.6 x 3.0 x 0.6mm



Posterior SG11/12



Posterior SG13/14



Posterior SG15/16

另備有Hu-Friedy相關器械目錄，歡迎來電洽詢